**MOPEDIA®**

MA R1830\_R1835 02 A\_ITA\_07 2020

**METTI IN PIEDI KOMPASS UP**  

---

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ITALIANO

## INDICE

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CE I Class Medical Device

REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL  
of 5 April 2017 concerning medical devices

## 1. CODES

### KOMPASS UP -STAND UP LIFT

RI830 - LINAK actuator max load 180 Kg

RI835 - TiMOTION actuator max load 180 Kg

## 2. INTRODUCTION

Thank you for choosing a patient stand up lift KOMPASS UP by MOPEDIA range, manufactured by MORETTI S.p.A. Their design and quality are a guarantee of comfort, safety and reliability. KOMPASS UP by MOPEDIA stand up lifts by MORETTI S.p.A. are designed and built to meet all your demands for a practical, correct use. This user manual provides for some suggestions as to how correctly use the lift you have chosen and offers a lots of valuable advices for your safety. Please read carefully through the manual before using the patient lift/stand up. Should you have any queries, please contact your retailer for advice and assistance.

## 3. INTENDED USE

Patient stand up lift KOMPASS UP by MOPEDIA is intended to stand up and move people with motor difficulties. An assistant operates the machine. The stand up lift cannot be used by one person only.

If you should have any doubts please contact your dealer who will help and advice you correctly.

### WARNING!

- Do not use the product for a purpose not indicated in this manual.
- Moretti S.P.A. declines all responsibilities for any consequences resulting from an incorrect use of this product and from unauthorized alteration to the frame of the product.
- The manufacturer reserves the right to change the information contained in this document without previous notice.

## 4. DECLARATION OF CONFORMITY

Moretti S.p.A. declares under its sole responsibility that the product made and traded by Moretti S.p.A. and belonging to the group of PATIENT LIFTS - MOPEDIA complies with the provisions of the regulation 2017/745 on MEDICAL DEVICES of 5 April 2017.

For this purpose, Moretti S.p.A. guarantees and declares under its sole responsibility what follows:

1. The devices satisfy the requirements of general safety and performance requested by the Annex I of regulation 2017/745 as laid down by the Annex IV of the above mentioned regulation.
2. The devices ARE NOT MEASURING INSTRUMENTS.
3. The devices ARE NOT MADE FOR CLINICAL TESTS.
4. The devices are packed in NON-STERILE BOX.
5. The devices belong to class I in accordance with the provisions of Annex VIII of the above mentioned regulation
6. Moretti S.p.A. provides to the Competent Authorities the technical documentation to prove the conformity to the 2017/745 regulation, for at least 10 years from the last lot production.

**Note:** Complete product codes, the manufacturer registration code (SRN), the UDI-DI code and any references to used regulations are included in the EU declaration of conformity that Moretti S.p.A. releases and makes available through its channels.

To guarantee greater safety, Moretti SpA followed the EN ISO 10535:2006 guidelines relating to lifters.

## 5. GENERAL WARNINGS

- For proper use of the device refer to this manual.
- For proper use of the device consult your doctor or therapist.
- Keep the packaged product away from any source of heat because the packaging is made of cardboard.
- The useful life of the accessory is determined by wear of parts that cannot be repaired and/or replaced.
- Be careful when children are nearby.
- The user and/or the patient will have to report any serious accident that have occurred related the device to the manufacturer and appropriate authority of the State which the user and/or patient belongs to.

### WARNING!

- Always pay close attention to the presence of moving parts that could cause limb entrapment and injury.

## 6. SYMBOLS

Product code

Unique Device Identification

CE mark

Manufacturer

Batch Lot

Read the instruction manual

Medical Device

Conditions of disposal

Attention, should read the instructions

Production dates

Waste Electrical & Electronic Equipment (WEEE)

Protection against dust and water splashes

B Type

## 7. GENERAL DESCRIPTION


The electric stand up lift KOMPASS UP helps people with motor difficulties to stand up or to be transferred. Ideal for wheelchair users, to stand up. It is the perfect aid to transfer and to stand up the patient from a sitting position with maximum load of 180 kg.

Perfect for those people who need assistance to stand up. Thanks to the width of the opening base the professional work is much easier and the transfers safer. The feet support is fixed and covered by a non-slip felt. The support for the legs is height adjustable. The control unit is provided with security systems.

Thanks to its control unit the electric stand up lift KOMPASS UP guarantees a comfortable and safe verticalization: the control unit monitors the battery charge level and the actuator conditions.

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




## 9. ASSEMBLY INSTRUCTIONS




### 9.1 Frame assembly R1830- R1835

1. Unscrew the screw at the base of the mast. (Fig.1, step 1)
2. Fit the mast into the provided cavity of the base. (Fig.1, step 2)
3. Fit the screw that blocks the mast into the base and screw it well. (Fig.1, step 3)
4. Connect the boom to the mast with the provided screws.(Fig.2, step 1-4)
5. Fit the handle into the mast and block it with the 2 screws. (Fig.3, step 1)
6. Fit the leg support to the mast and block it with the 2 screws. (Fig.3, step 2)
7. Fit the lever for opening the base and block it with its screw. (Fig.3, step 3)
8. The structure will be completely assembled. Install the actuator suitable for the chosen patient lift.

### 9.2 Assembly LINAK electrical kit (for R1830)

1. Secure the control unit support bracket to the mast with the provided screws, using the fitted cavities. (Fig.4)\*
2. Insert the battery in the appropriate place.
3. Fit the elevation actuator between the boom and the mast and secure it firmly with screws and bolts. (Fig.5)\*
4. Connect the electrical connector of the actuator to the control unit socket "1" - make sure it is fully inserted. Connect the remote control electrical connector to the large socket on the control unit - make sure the guide lug seats properly and then push it all the way in. (Fig.6)\*


\*images for illustrative purpose only

### 9.3 Assembly TIMOTION electrical kit (for R1835)

1. Secure the support bracket of the control unit to the mast with the provided screws, using the fitted cavities. (fig.7)\*
2. Fit the control unit on the bracket, letting it slide from top to bottom and blocking it with the screws. (Fig.8)\*
3. Fit the battery into the control unit. (Fig.9)\*
4. Remove the two actuator mounting bolts from the mast boom. Fit the actuator into it and secure it with the two bolts previously removed. Cover the bolt heads with the provided plastic cap. (Fig.10)\*

5. Connect the actuator connector to the central socket of the control unit and the remote control to the right socket. (Fig.11-12)\*

\*images for illustrative purpose only

## 10. BEFORE USE

KOMPASS by MOPEDIA stand up lifts are intended to move people with motor difficulties for short transfers as from a bed to a wheelchair or a commode chair and vice versa. The stand up lift cannot be used to transport the patient.

- Before lifting a patient make sure the base legs are opened to the maximum width possible, that the opening lever is secured in the opened position and the rear wheels brakes are blocked.
- While moving the patient, keep an eye out for obstacles that might cause the lift to tip over.

⚠ Always check that the sling is in good condition. Do not use slings showing signs of damage or wear

## 11. INSTRUCTIONS FOR USE

BEFORE USING THE ELECTRICAL LIFT YOU ARE RECOMMENDED TO CHARGE THE BATTERY FOR ABOUT 6/8 HOURS. AFTER THAT MAKE SURE THE BATTERY IS CHARGED BEFORE EVERY USE.

Before lifting a patient, open the base legs to the maximum width possible, gently pulling the opening lever towards you and secure it in the open position

1. Block the rear wheel brakes.
2. Keep the patient centered between the base legs and always turned towards the assistant.
3. Adjust the leg support position in height and width using the regulation knobs.
4. Put the sling on the patient, adjusting the abdominal band and the underarm belts.
5. Hook carefully the sling's straps to the lifting beam, see picture (The central hook is to be used with the pelvic sling. The middle hooks are to be used by persons < 1.70m tall. The hooks at the ends of the arm are to be used by persons > 1.70m tall)

**NOTE:** A pelvic sling is recommended for patients with a sturdy build or low mobility. Thanks to the pelvic belts the underarm pressure decreases and the verticalization becomes more comfortable for the patient

**NOTE:** If necessary, for tall persons you can hang the sling on the arm middle hooks, passing the belts through the ends hooks. In this case, the underarm belts become shorter and let the patient's chest approach the arm and incline toward it.

6. Press the up-facing arrow key to lift the patient.
7. ALWAYS MONITOR THE PATIENT'S POSITION THROUGHOUT THE LIFTING PROCESS!
8. Support the patient until the desired height, keeping him facing his assistant.
9. Release the rear wheels and move the patient.
10. Before lowering the patient, block the rear wheels again.
11. To lower the patient, press the down-facing arrow key.
12. ALWAYS MONITOR THE PATIENT'S POSITION THROUGHOUT THE LIFTING PROCESS!

**WARNING!**


- Always help the patient during the lowering process so as to maintain the correct seated position!
- While lowering the patient, please remove any barrier that can impede the downward movement. If the case, lift up the patient again and remove what causes the impediment for lowering.

**Automatic anti-crushing system**

The electrical kit is equipped with an automatic anti-crushing system, which stops the lowering of the lifting boom if there is an accidental contact between the patient and any barrier.

**Emergency stop**

The control unit is equipped with a red "STOP" button which instantly deactivates all controls when pressed. To restore normal operation, turn the button clockwise and pull it outwards.

#### Mechanical lowering of the actuator

If using the manual/mechanical lowering, please make sure that it will be done in total safety conditions for the patient. In order to activate the manual/mechanical lowering, pull downwards the red ring placed in the upper part of the actuator and turn it counterclockwise. Once finished the lowering process, take the red ring back to its resting position.

#### Emergency keys

**RI830:** Underneath the red STOP button of the control box, there is a small button that can be used to low patients when the remote control fails to operate properly. This button can only be operated with the aid of a narrow tipped object like the tip of a pen. (Fig.13)

**RI835:** On the control unit there are two keys (one with an up-facing arrow and the other one with a down-facing arrow) that can be pressed instead of those ones on the remote control to complete all the operations when the remote control fails to work properly.

#### Indication of the battery state of charge (only RI835)

There is a 3-Led indicator for the battery state of charge. This indicator is placed on the remote control. 3 Leds on indicate a complete charge of the battery, whereas only one Led on indicates the battery is out of power and needs to be charged.

#### Recharging batteries

- Connect the adapter to the main socket and then to the socket on the control unit, after removing the protection cap.
- On the charger adapter there is a led that during the charging process becomes orange color. Charge the battery until the led becomes green color (RI835 only)
- Before activating the stand up lift again, make sure to disconnect the battery charger!
- Charge the batteries for at least 6 hours before using the stand up lift.

## 12. MAINTENANCE

The MOPEDIA devices by MORETTI are checked carefully and supplied with EC mark, once launched on the market. For patient and physician's safety we recommend to let check the device by the manufacturer or by laboratory authorized every 2 years. In case of repair, use only original spare parts and accessories.

- It is extremely important to check the lift before each use.
- More specifically, please check that all screws and bolts are well tightened.
- Please, check the lift can move freely with the wheels unlocked by the brakes.
- Check the rear wheel brakes function correctly.
- It is very important to regularly check the connection pins and stop rings are correctly located (the pins connect the actuator to the upper and lower part of the lift, and the lifting beam where you hook the cradle).
- The mechanical parts must be lubricated with a drop of oil in order to facilitate their sliding and resistance to wear. In particular, pay attention to the mechanical retractor for the legs of the patient lift. Lubricate as indicated above every two or three months.

## 13. CLEANING AND DISINFECTION

### 13.1 Cleaning

To clean the cradle, handle and the structure use a solution of alcohol and neutral soap (in equal parts) mixed with water.

### 13.2 Disinfection

If you need to disinfect the product use a common disinfectant neutral detergent.

**Note:** Never use acids, alkalines or solvents such as acetone or diluent. Keep away from light, moisture, heat.

## 14. CONDITIONS OF DISPOSAL

In case of disposal don't use the inserting container for municipal waste. We recommend to dispose the kit in the appropriate disposal areas for recycling.

### 14.1 Warnings for the proper disposal of the product, as required by the European Directive 2012/19/UE

At the end of its service life, the product must not be disposed of together with normal municipal waste. It must be delivered to the municipal sorted waste collection centre or to the resellers, in case they provide for this service. Sorted waste disposal reduces the environmental and health consequences of an improper disposal. Furthermore, the materials which make up the product may be recycled, thus resulting in energy and resource savings. In order to highlight the compulsory sorted disposal of the electro-medical equipment, the product itself bears the barred bin symbol.

### 14.2 Treatment of dead batteries (directive 2006/66/CE)

The symbol on the product indicates that the batteries must not be treated like common household waste. Make sure the batteries are properly disposed; this prevents the possible environmental and health consequences of an improper disposal. Recycling helps to preserve natural resources. Deliver dead batteries to a proper recycling centre. For more detailed information regarding the disposal of dead batteries or of the product itself, please contact the municipality, the local waste disposal centre or the reseller from whom you purchased the device.

## 15. SPARE PARTS ACCESSORIES

RP893 Body harness with dorsal abdominal band with Velcro closure

RP896 Body harness with dorsal abdominal band with Velcro closure and pelvic belts

For spare parts and accessories, please always refer to the MORETTI master catalogue.

**ONLY USE ORIGINAL MORETTI ACCESSORIES. MORETTI S.P.A. COMPANY SHALL NOT BE LIABLE FOR ANY INCIDENT CAUSED BY IMPROPER USE.**

## 16. TECHNICAL FEATURES

### 16.1 Technical specifications

- Powder coated steel frame
- Platform covered with felt
- Knee pad adjustable in height and width
- Castors with brake and rotational block
- Leg opening lever with safety lock
- Comfortable handle for the healthcare professional
- Hooks designed specifically to hold the sling in complete safety
- Autonomy: 100 to 200 lifts with fully charged battery
- Total weight: 40 Kg
- Maximum loading weight: 180 Kg

### 16.2 Frame dimensions - RI830-RI835

#### Legend

A Sagittal dimension at max.verticalization: 29 cm

B Arm length: 81 cm

C Maximum stroke: 80 cm

D Minimum arm height: 100 cm

E Maximum arm height: 180 cm

F Maximum base length: 110 cm

G Distance between front castors and mast: 70 cm

H Minimum internal width: 44 cm



I Minimum external width: 59 cm

L Maximum internal width: 87 cm

M Maximum external width: 100 cm

O Castors height: 10,5 cm

P Legs height: 11,5 cm

## 17. WARRANTY

Every Moretti product is guaranteed for 2 (two) YEARS from the purchase date, except the restrictions specified below. The warranty shall not apply in case of damage caused by improper use, abuse, alterations of the product or any other use not in compliance with the user manual. The proper use of the product is indicated in the user manual. Moretti S.p.A. shall not be liable for any damages or personal injuries due to incorrect assembly/improper use of the device, not compliant with the user manual. Moretti does not guarantee its products against damages or faults in the following cases: natural disasters, non authorized maintenance or reparations, damages caused by electrical supply problems (if provided), use of parts or components not provided by Moretti, failure in accomplishing the guidelines and instructions for use, non authorized alterations, damages due to the transport (different from Moretti original delivery) or failure in maintaining the device as indicated in the manual. The components subject to wear are not covered by this warranty if the damage is due to the normal use of the product.

### 17.1 Warranty for rechargeable battery (if provided)

The original batteries and the spare batteries are covered by a (ninety) 90-days warranty as far as their performance is concerned and by a (six) 6-months warranty as far as manufacturing faults are concerned or as required by law. Warranty on fully charged batteries left in an unused state for more than (three) 3 consecutive months is automatically void. Warranty on fully discharged batteries left in an unused state for more than (three) 3 consecutive days is automatically void.

### 17.2 Repairing Warranty repair

If a Moretti product shows material or manufacturing faults during the warranty period, Moretti S.p.A., together with the client, will consider whether the product fault is covered by warranty. Moretti S.p.A., at its unquestionable discretion, may replace or repair the item under warranty in a Moretti reseller point or at its own qualified offices. Labour costs for repairing the product may be at Moretti's expense if it is determined the repair is under warranty. A repair or replacement does not renew nor extend the warranty.

### Repair not covered by warranty

A product, not covered by warranty, may be returned for reparation only after receiving prior authorization from Moretti client service. Labour and delivery costs for an out-of-warranty reparation will be fully at the client's or reseller's expense. Reparations of out-of-warranty products are guaranteed for 6 (six) months, from the day the client receives the repaired product.

### Non defective device

The client is informed if, after examining and testing a returning product, Moretti decides the product is not defective. The product will be sent back to the client and all the return delivery costs will be at his expense.

### 17.3 Spare parts

The original Moretti spare parts are guaranteed for 6 (six) months from the day you receive them.

### 17.4 Exempt clauses

Except what clearly specified in this warranty and in accordance with the law, Moretti S.p.A. does not offer its clients any other declaration, warranty or condition, expressed or implicit, including any declaration, warranty or condition of merchantability, of non violation or non interference, of suitability for a specific aim. Moretti S.p.A. does not guarantee that the use of Moretti products will be continuous and faultless. The duration of possible implicit warranties that may be applied in accordance with the law is limited to the warranty period.

Some states or countries do not allow limitations to the duration of an implicit warranty, or the exoneration, or the limitations of accidental or indirect damages related to products for consumers. In such states or countries some exonerations or limitation of this warranty may not be applied to the user. This warranty is subject to modifications without any prior notice

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